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Table of contents

▪ Used symbols.....	17
▪ Your LED illumination system consists of the following components.....	18
▪ Device setup.....	20
1. Charging the battery pack.....	20
2. Control elements on the unit.....	20
▪ Installation of the yellow filter.....	21
▪ Connecting the control unit with the cable/spotlight	21
▪ Attaching the spotlight to an Orascope loupe	22
▪ Removing the belt clip	22
▪ Cleaning and disinfection of the spotlight.....	23
▪ Disposal instructions.....	24
▪ Troubleshooting.....	24
▪ Technical data	25
▪ Safety and warning instructions	26
▪ Limited manufacturer's warranty	27

The LED illumination systems " **Ultra Smart PLUS ZWELL** " and "**Ultra Smart PLUS STUDY**" are portable medical devices. They are used as examination lights by dental professionals (dentists, dental hygienists, dental assistants) to illuminate the oral cavity. This illumination systems will be applicable in all kinds of examination rooms in which dental work could be carried out, also very often in combination with a loupe.

Used symbols

Symbols	Function / Meaning
	Medical Device – indicates that the item is a medical device
	Order number / Item number
	Serial number
	Manufacturer
	Year of manufacture
	Product is compliant with EU - Medical Devices Directive/Medical Devices Regulation
	Not for the general waste
	Caution (Check accompanying documents)
	Do not look directly into the light
	Consult instructions for use
	Protection class II device


Thank you for choosing our LED headlights. This instruction manual is intended to help you to know better the range of functions of your new LED headlight and the optimal use of it.

Please read this instruction manual fully before using the LED headlight for the first time! Please keep this instruction manual for future use.


Your LED illumination system consists of the following components:

1. Spotlight with magnetic system
2. Magnetic connector cable
3. Control unit with belt clip
4. Power supply / Battery charger (please follow the manufacturer's operating instructions that are also enclosed)
5. Ultra Smart filter kit / content:
 - a) Filter clamp
 - b) Yellow filter (prevents premature curing of composite materials)
 - c) Transparent glass as protection against contamination
6. TTL adapter (incl. 2 pcs. Phillips screws)
7. Cable guide (to attach on the temple arm)
8. Cleaning brush

1. Spotlight with magnetic system

2. Magnetic connector cable

3. Control unit with belt clip

4. Power supply / Battery charger

5./6./7./8. Accessories

Yellow filter
incl. filter clamp

Transparent
glass

TTL adapter
incl. 2 pcs. Phillips screws
(Illustration similar)

no illustration

Cable guides

Cleaning brush

Device setup

1. Charging the battery pack

The lithium-ion battery pack in the control unit is not fully charged upon delivery. Prior to first use, please proceed as follows:

- Attach the battery charger to the control unit.
- Attach the charger to the power outlet.
- The indicator of the charger shines "red" when charging.
- Once the battery pack is fully charged, the charger indicator light turns to "green" and switches automatically to trickle charging. The control unit must be separated from the charger and is ready for use.

IMPORTANT! It will need three to four times of using and recharging of the control unit to reach the maximum capacity of the battery pack.

2. Control elements on the unit

Ultra Smart PLUS zwei	Ultra Smart PLUS STUDY
a) Turning knob = Brightness control Controls the light intensity of the spotlight	a) Rocker Switch = Brightness control Controls the light intensity of the spotlight
Level 1 = lowest intensity Level 2 = medium intensity Level 3 = highest intensity	Level 0 = device is off Level 1 = lowest intensity Level 2 = highest intensity
Stage „off“= The spotlight is off but the device is at "monitoring mode", i.e. the present capacity of the accumulator is being checked regularly. An acoustic signal indicates the dropping of the capacity below the minimum level.	

(Light's duration per stage and light's intensity → refer to the technical data)

Ultra Smart PLUS zwei	Ultra Smart PLUS STUDY
b) Battery indicator	b) Battery indicator
If the remaining capacity of the battery drops to the minimum level (about 40 %) the indicator light at the control unit flashes "red". When it drops further down to 30 % the flashing light turns into a constant light. Additionally, an acoustic signal sounds every 10 minutes.	Illumination ring = shines (when switched on) permanently green when charged. If the capacity drops below 30 % of the total capacity, the ring changes to permanent red. Approximately 10 minutes before complete discharge, the ring then flashes red.

Charge the accumulator as described in "Charging the battery pack".

IMPORTANT! The battery capacity depends always on the chosen light intensity!

Installation of the yellow filter

1. Loosen the hinge screw, slide the yellow filter onto the spotlight with the hinge up and screw it tight.
2. The filter is ready to use and could be fold down to the spotlight if needed.

Connecting the control unit with cable/spotlight

1. Connect the plug of the magnetic connector cable with the control unit.

IMPORTANT! The red dots at the plug as well as at the socket must be aligned accurately!

2. The magnetic connector cable as well as the spotlight contain a magnet for connecting the cable to the spotlight. **The magnets are attracted automatically so do not force the magnetic connector cable onto the spotlight!**

3. In case of separating the magnetic connector cable from the control unit never unplug by pulling direct at the cable but only at the ribbed section of the plug instead!

IMPORTANT! Do not overstretch the magnetic connector cable, e.g., do not wrap it around the control unit.

Attaching the spotlight to an Orascopeptic loupe

1. To attach the spotlight the TTL adapter must be mounted at the bridge of the frame with the two provided Phillips screws.

2. The spotlight (screw down) is then inserted into the TTL adapter device.
3. To remove it from the TTL adapter device (consists of a high quality and robust material) just pull the spotlight vertically up.

Removing the belt clip (only Ultra Smart PLUS ZWELL)

1. Pull the release locking mechanism downwards and slide the unit upwards

Cleaning and disinfection of the spotlight

1. The spotlight comprises the highest quality components and therefore it must be handled with the utmost care!
2. Magnet points must be kept clean at all times with the enclosed brush and a (disinfectant) cloth that has been moistened with a suitable disinfectant solution. Wipe dry the cleaned contacts. Any contamination could affect the position of the spotlight or release a light interruption.
3. In order not to scratch the lens surface impurities should be removed by using a soft, clean, and dry cloth only!
4. **Important:** The spotlight is not waterproof! Do not immerse the light in water, disinfectant liquids or ultrasonic cleaning liquids!
5. No autoclaves or chemiclaves may be used for sterilization.
6. The use of glutaraldehyde, methanoic acid (formic acid) and agents containing iodine may cause damages to the spotlight!
7. Please use a damp cloth for disinfection only (such as „MinutenWipes“/ „PlastiSept-Wipes“ by Alpro).

In the case of utilizing other moist disinfecting cloths, make sure that your moist disinfecting cloths:

- a) are suitable for the surfaces of medical products that are sensitive to alcohol (e.g., dental chairs, dental instrument tables, instrument surfaces, etc.),
- b) do not contain any form of 1-propanol, aldehyde, phenol, formic acid, or any quaternary ammonium compound.
- c) are bactericide, fungicide and virucide

before using them for the first time.

Disposal instructions

In accordance with the legal requirements (ElektroG: Waste electrical equipment) we take back our equipment. It will be either recycled or disposed by a recycling company in accordance with legal requirements.

Troubleshooting

Problem	Solution
The control unit is turned on and attaches with the magnetic connector cable, but the spotlight is off.	▪ Check whether there is any contamination (dirt or metal particle) on one or both magnets and clean if necessary as described on page 23. ▪ Make sure that the light is properly connected to the control unit and the magnetic connector cable. ▪ Recharge the control unit.
The light is off / light has low luminosity	▪ Upgrade the light intensity regulator. ▪ The battery pack might be too weak → Recharge the control unit.
Light interruption	▪ Check whether there is any contamination (dirt or metal particle) on one or both magnets and clean if necessary as described on page 23.
Power supply to the battery charger is interrupted before the end of the load cycle (for example power failure).	▪ As soon as the power supply is re-established, the charge cycle will restart.
Indicator light of battery charger remains "red" (even after more than 8 hours of charging) and does not switch to "green".	▪ Return illumination system including battery charger for inspection to Sigma Dental.
Indicator light of battery charger remains "red" even after being reconnected from the control unit.	▪ Briefly disconnect battery charger from the power grid and return the equipment for inspection to Sigma Dental.

Technical data

	Ultra Smart PLUS ZWEI	Ultra Smart PLUS STUDY
Weight of spotlight	5,6 g	5,6 g
Weight of control unit	170 g	82 g
Diameter of spotlight	20,7 mm	20,7 mm
Dimensions of control unit	62 mm x 97 mm x 26 mm	56 mm x 82 mm x 24 mm
Voltage / Charge capacity	3,7 V DC / 6400 mAh	3,7 V DC / 2400 mAh
Battery type	Li Ion	Li Ion
Battery charging time	approx. 7 hours	approx. 3,5 hours
Luminous period	Intensity Level 3 (max.) approx. 13 hrs Intensity Level 2 (med.) approx. 22 hrs Intensity Level 1 (min.) approx. 32 hrs	Intensity Level 2 (max.) approx. 10 hrs Intensity Level 1 (min.) approx. 15 hrs
Efficacy	approx. 80.000 LUX measured at a distance of 30 cm	approx. 60.000 Lux measured at a distance of 30 cm
Battery charger type	Mascot 2241000074	Mascot 2241000073
Power output of battery charger	Max. 4,2 V / 1,3 A	Max. 4,2 V / 1,3 A
Power supply voltage of battery charger	100 – 240 V / 50 – 60 Hz / max. 0,2 A	100 – 240 V / 50 – 60 Hz / max. 0,2 A
Ambient air temperature for operating and storage	0 to 40° C	0 to 40° C
Interference resistance information	Non-life-sustaining equipment (IEC 60601-1-2)	Non-life-sustaining equipment (IEC 60601-1-2)
Photo-biological safety	DIN EN 62471 (IEC/TR 62471) Risk class 2	DIN EN 62471 (IEC/TR 62471) Risk class 2

Safety and warning instructions

1. To avoid the risk of explosion, use the control unit with the supplied charger Mascot Li-Ion 2241.
2. To prevent the risk of permanent injury to health (such as blindness) do not focus the light on eyes injuries must wear safety glasses.
3. The control unit must be used only for the intended purpose.
4. By no means the housing of the battery recharger or the control unit should be opened. Do not disassemble or modify the equipment at all. The improper use of the control unit may cause excessive heating, explosion or igniting of the battery pack and thus may lead to serious injuries. Please make sure that the safety instructions below are followed:
 - Do not open
 - Do not modify the equipment
 - Do not expose to fire or excessive heat (above 70° C)
 - Do not use outdoors
 - Do not short-circuit (do not connect contacts via a conductor or reverse polarity)
 - Protect against fluids or humidity
 - Do not use damaged equipment
 - Follow temperature range from 0° - 40° C for storage and operating.
5. Disconnect the battery charger from the control unit when not in use.
6. Disconnect the battery charger from the wall socket when not in use.
7. Never use any faulty battery chargers or control units.
8. The battery charger must not be used outdoors.
9. Do not use the battery charger or control unit if the housing or the plug are damaged.
10. In case the battery charger or control unit become increasingly hot, have a burnt smell, change colour or shape – please switch off the battery charger immediately and disconnect it from the power grid. By no means the equipment may be used any longer. In this case it must be returned to Sigma Dental directly.
11. In case acid leaks from the control unit and come into contact with skin or clothing, rinse immediately under running water to prevent skin irritations.
12. In case of a serious incident which occurs in connection with the device, manufacturer and competent authority (member state of user) must be informed.
13. Do not keep within reach of children.

14. Temperature variations may cause of condensate. To avoid malfunction the device should be adjusted to room temperature for at least 2 hours before use.

15. **CAUTION:** Do not look directly into the light. Looking directly into it from a short distance could result in severe damages for the eyes.

LED safety class according to IEC 62471: risk class 2

Limited manufacturer's warranty

This limited manufacturer's warranty by Sigma Dental Systems – Emasdi GmbH („Sigma Dental“) shall be applicable to all European countries unless a separate guarantee is in existence locally. Pursuant to the following contract terms and conditions, Sigma Dental shall warrant that, at the time of the original purchase, the equipment was free of defects in terms of material, construction and implementation:

1. This limited warranty shall be applicable to the end purchaser of the equipment („Customer“). Any statutory claims of customer or of customer towards seller / dealer shall neither be excluded nor restricted by same.
2. The warranty period shall be twelve (12) months following the purchase of the equipment by the original customer. In case the equipment is resold or otherwise changes owner / user, the warranty shall be valid for the remainder of the 12-month period. These terms and conditions shall not be affected by such change. This limited manufacturer's warranty shall be valid and enforceable only in countries of the European Union.
3. During the warranty period, any defective equipment (at the sole discretion of Sigma Dental) shall be repaired or replaced by Sigma Dental or an authorized customer service. The customer shall receive either the repaired equipment or a fully functional replacement item. Any replaced parts or accessories shall become the property of Sigma Dental.
4. There shall be no extended or renewed warranty period for repaired or replaced equipment.
5. This limited warranty shall not be applicable to any equipment painted or otherwise personalized.
6. Any consumables (such as battery pack) shall have a warranty period of 6 months after first purchase.

7. This limited warranty shall not be applicable to normal symptoms of wear and tear. Furthermore, this limited warranty shall not be applicable if:
- a) the defect was caused by non-compliance with the operating instructions, improper handling, wetness, moisture, extreme heat of climate conditions or short-term fluctuations of such effects, corrosion, oxidation, unauthorized access or connection attempts, unauthorized opening and repairs, repair attempts with unauthorized spare parts, operator error, improper installation, accidents, natural disasters, spillage of beverages or food, chemical effects or other external effects outside the realm of influence of Sigma Dental (including defects in consumables, such as battery packs, which necessarily have a limited lifespan), unless the defect is directly caused by an error in material, construction, or manufacture;
 - b) the customer does not report the defect to Sigma Dental or the authorized customer service within thirty (30) days after the defect occurred;
 - c) the serial number of the equipment has been removed, scratched out, crossed out or modified or the serial number is illegible;
 - d) the defect was caused by the use of the equipment in combination with a part not authorized by Sigma Dental, or by connecting the equipment to such a part, or by using the equipment in some other way in violation of its intended use;
 - e) the defect was caused by short-circuiting the battery pack, or by breaking the wrapping of the battery pack, or manipulation same, or by using the battery pack for a device other than the intended device.
8. This limited warranty shall represent the only and exclusive legal remedy of the customer against Sigma Dental, and the only and exclusive liability of Sigma Dental toward the customer for defects or functional problems of the equipment. This limited warranty shall replace all other guarantees and warranty statements, whether provided orally or in writing, and whether based on (not necessarily forcefully applicable) statutory regulations, contract provisions or unauthorized acts. Under no circumstances Sigma Dental shall be liable for any ancillary, subsequent or direct damages, costs or expenditures. If the customer is a legal entity, Sigma Dental shall not be liable for any indirect damages, costs or expenditures.
9. Upon assertion of any warranty claim, the original receipt, including the date of sale, must be enclosed.
10. Any changes to this limited warranty shall require the prior written approval by Sigma Dental.

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