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Table of contents

▪ Used symbols.....	17
▪ Components of LED illumination system.....	18
▪ Device setup.....	19
▪ Installation of the yellow filter.....	20
▪ Attaching the spotlight to an Orascoptic loupe.....	20
▪ Connecting the spotlight with the control unit.....	21
▪ Activation and deactivation of the spotlight / reduction of the light intensity.....	21
▪ Cleaning and disinfection of the spotlight.....	23
▪ Disposal instructions.....	23
▪ Troubleshooting.....	24
▪ Technical data.....	25
▪ Safety and warning instructions.....	26
▪ Limited manufacturer's warranty.....	28

The LED illumination system "Ultra Smart^{plus} CU•BE touch" is a portable medical device. It is used as an examination light by dental professionals (dentists, dental hygienists, dental assistants) to illuminate the oral cavity. This illumination system will be applicable in all kind of examination rooms in which dental work could be carried out, also very often in combination with a loupe.

Used symbols

Symbols	Function / Meaning
	Medical Device – indicates that the item is a medical device
	Order number / Item number
	Serial number
	Manufacturer
	Date of manufacture
	Product is compliant with EU Medical Devices Directive/Medical Devices Regulation
	Not for the general waste
	Caution (Check accompanying documents)
	Do not look directly into the light
	Consult instructions for use
	Protection class II device

Many thanks for choosing the wireless "Ultra Smart^{plus} CU•BE touch" illumination system. This instruction manual is intended to help you to know the range of functions of your new illumination system and the optimal use of it.

Please read this instruction manual carefully before using the illumination system for the first time and keep this instructions for future reference.

Your wireless LED illumination system consists of the following components:

1. Spotlight with magnet
2. Two control units with rechargeable battery pack, magnet and touch function
3. USB charging stick with magnet
4. USB power plug
5. Ultra Smart filter kit / content:
 - a) Filter clamp
 - b) Yellow filter (prevents premature curing of composite materials)
 - c) Transparent glass as protection against contamination
6. TTL adapter (incl. 2 pcs. Phillips screws)
7. Cleaning brush

1. Spotlight with magnet

2. Two Control units with rechargeable battery pack, magnet and touch function

3. USB charging stick with magnet

4. USB power plug

- 5./6./7. Accessories

Device setup

Charging of the control unit


The battery pack in the control unit is not fully charged upon delivery. Prior to first use, please proceed as follows:

- Connect the USB power plug with the power outlet.
- Then connect the USB charging stick with the power plug (indicator lamp "green"). USB charging stick also connects to any USB port 5 V/200 mA.
- Thereafter connect the control unit and the USB charging stick to start charging.
- The indicator of the USB charging stick shines "red" while charging.
- Once the battery pack is fully charged (or ready for use), the indicator light of the USB charging stick turns to "green".
- The control unit will be separated from the USB charging stick and it can be used now.

IMPORTANT!

- The final capacity of the battery pack is not reached before the equipment has been used and discharged three to four times, and then recharged again.
- If the control unit is not used for a long period, a new charging is required.
- Magnet points must be kept clean at all times with the enclosed brush and a (disinfectant) cloth that has been moistened with a suitable disinfectant solution (see page 23). Wipe dry the cleaned contacts. Any contamination could affect the position of the light or release a light interruption.

Installation of the yellow filter

1. Loosen the screw, slide the yellow filter onto the spotlight and screw it tight.
2. The filter can be fold down now in front of the spotlight if necessary.

Attaching the spotlight to an Orascopeptic loupe

1. To attach the spotlight to the Orascopeptic loupe, the TTL adapter is mounted over the bridge on the loupe frame with the two supplied Phillips screws.

2. The spotlight is then inserted into the TTL adapter device. Hold the spotlight over the adapter so that the screw of the spotlight points downwards. Now insert the screw into the hole of the adapter.

Connecting the spotlight with the control unit

1. The control unit as well as the spotlight include a magnet.
2. If the spotlight is fixed in the TTL adapter, adapt the fully charged control unit with the magnet point to the magnet point of the spotlight. The resulting magnetic force pulls the control unit towards the spotlight.

3. As soon as the control unit is attached, the spotlight starts to shine (at the highest light intensity – level 2) and can be used immediately.
4. At the end of the process, disconnect the control unit from the spotlight.
5. To avoid possible damage, please do not put the two components together with additional pressure.

Activation and deactivation of the spotlight / reduction of the light intensity

The control unit has two touch-sensitive surfaces, marked with the word "touch" on top and the serial number on the bottom. After connecting the control unit and spotlight (the spotlight automatically starts at the highest light intensity – level 2), please wait about three seconds. The control unit adjusts itself during this time. The touch function is then activated.

1. Reduction of the light intensity (to level 1): Tap once briefly on one of the touch-sensitive surfaces of the control unit.
2. Switching off the spotlight: briefly tap again one of the touch-sensitive surfaces.
3. Switching the spotlight on again: Tap again one of the touch-sensitive surfaces. The spotlight automatically start at light level 2.

4. At the end of the process disconnect the control unit from the spotlight

IMPORTANT!

When installing the TTL adapter, control unit and spotlight, be sure that the positioning is correct to avoid the control unit tilting or dropping off.

correct:

incorrect:

incorrect:

Cleaning and disinfection of the spotlight

1. The spotlight comprises the highest quality components and therefore it must be handled with the utmost care!
2. Magnet points must be kept clean at all times with the enclosed brush and a (disinfectant) cloth that has been moistened with a suitable disinfectant solution. Wipe dry the cleaned contacts. Any contamination could affect the position of the light or release a light interruption.
3. In order not to scratch the lens surface impurities should be removed by using a soft, clean and dry cloth only!
4. The spotlight is not waterproof! Do not immerse the light in water, disinfectant liquids or ultrasonic cleaning liquids.
5. No autoclaves or chemiclaves may be used for sterilization. The use of glutaraldehyde and agents containing iodine may cause damages to the spotlight!
6. Please use a damp cloth for disinfection only (such as „MinutenWipes“ / „PlastiSept-Wipes“ by Alpro).

In the case of utilizing other moist disinfecting cloths, be absolutely sure that your moist disinfecting cloths

- a) are suitable for the surfaces of medical products that are sensitive to alcohol (e.g., treatment units, dental instrument tables, instrument surfaces, etc.),
- b) **do not** contain any form of 1-propanol, aldehyde, phenol, any quaternary ammonium compound (such as cavidol) or formic acid
- c) are bactericide, fungicide and virucide.

before using them for the first time.

Disposal instructions

In accordance with the legal requirements (ElektroG: Waste electrical equipment) we take back our equipment. It will be either recycled or disposed by a recycling company in accordance with legal requirements.

Troubleshooting

Problem	Solution
The control unit is connected but the spotlight is off.	▪ Check whether there is any contamination (dirt or metal particle) on one or both magnets and clean if necessary as described on page 23. ▪ The resulting magnetic force pulls the control unit towards the spotlight. Do not press with force the spotlight into the magnet of the control unit. ▪ Make sure, that the spotlight is properly connected to the control unit. ▪ Recharge the control unit.
The spotlight is off / spotlight has low luminosity.	▪ Recharge the control unit.
Power supply to the USB charger is interrupted before the end of the charging cycle.	▪ As soon as the power supply is restored, the charging cycle will restart.
Indicator light of USB charging stick does not switch from "green" to "red", although the control unit is connected.	▪ Check whether there is any contamination (dirt or metal particle) on one or both magnets and clean if necessary as described on page 23. ▪ Make sure, that the control unit is properly connected to the USB charging stick. ▪ Battery pack in the control unit is already fully charged.
Indicator light of USB charging stick remains "red", although no control unit is connected	▪ Disconnect USB charging stick from USB power plug and connect it again.
Indicator light of USB charging stick does not light up "red" or "green".	▪ Check whether the USB power plug works. If necessary test with another USB port.
Control unit is getting warm.	▪ Remove control unit from the spotlight or USB charging stick and send it to Sigma Dental (address on the back).
Lighting period of the control unit is considerably reduced.	▪ Check the USB connection. ▪ Replace control unit, because battery pack capacity is decreasing.
Light interruption	▪ Check whether there is any contamination (dirt or metal particle) on one or both magnets and clean if necessary as described on page 23.

Technical data

Weight of spotlight	approx. 5,6 g
Weight of control unit	approx. 13,5 g
Diameter of spotlight	20,7 mm
Dimensions of control unit	31 mm x 23 mm x 18,5 mm (W x H x D)
Voltage / Charge capacity	3,7 V DC / 400 mAh
Battery type	Li Ion
Battery charging time	approx. 2 hours
Luminous period (step 2)	approx. 2 hours
Luminous period (step 1)	approx. 3 hours
Efficacy (step 2)	approx. 45.000 Lux
Efficacy (step 1)	approx. 30.000 Lux measured at a distance of 30 cm
Power output of USB charging stick	max. 4,2 V / 260 mA
Power output of USB power plug	max. 5 V / 500 mA
Input voltage USB power plug	100 – 240 V / 50 – 60 Hz
Ambient air temperature for operating and storage	0 to 40° C
Interference resistance information	Non life-sustaining equipment (IEC 60601-1-2)
Photo-biological safety	EN 62471 (IEC/TR 62471) Risk class 2

Safety and warning instructions

1. In order to avoid the risk of explosion, use the control unit with the supplied USB charger only.
2. To prevent the risk of permanent injury to health (such as blindness) do not focus the light on eyes. Persons at risk for eye injuries must wear safety glasses.
3. Battery pack in the control unit must be used only for the intended purpose (see page 16).
4. By no means the housing of the USB power plug, USB charging stick or control unit should be opened. Do not disassemble or modify the equipment at all. The improper use of the battery pack may cause excessive heating, explosion or igniting of the battery pack and thus may lead to serious injuries. Please make sure that the safety instructions below are followed:
 - Do not open
 - Do not modify the equipment
 - Do not expose to fire or excessive heat (above 50°C)
 - Do not use outdoors
 - Do not short-circuit (do not connect contacts via a conductor or reverse polarity)
 - Protect against fluids or humidity
 - Do not use damaged equipment
 - Follow temperature range from 0° – 40° C for storage and operating.
5. Disconnect the USB power plug from USB charging stick and also the control unit when not in use.
6. Disconnect the USB power plug from the wall socket when not in use.
7. Never use any faulty USB chargers or control units.
8. Do not use the USB power plug or control unit if the housing or the plug are damaged.
9. In case the USB power plug or control unit become increasingly hot, develop a smell of burning, change color or shape - please disconnect it from the power grid. By no means the equipment may be used any longer. In this case it has to be returned to Sigma Dental directly.
10. In case acid leaks from the control unit has come in contact with skin or clothing, rinse immediately under running water in order to prevent skin irritations.
11. In case of a serious incident which occurs in connection with the device, manufacturer and competent authority (member state of user) must be informed.
12. Do not keep within reach of children.

13. Temperature variations may cause of condensate. To avoid malfunction the device should be adjusted to room temperature for at least 2 hours before use.

14.

CAUTION: Do not look directly into the light. Looking directly into it from a short distance could result in severe damages for the eyes.

LED safety class according to IEC 62471:risk class 2

Limited manufacturer's warranty

This limited manufacturer's warranty by Sigma Dental Systems – Emasdi GmbH („Sigma Dental“) shall be applicable to all European countries unless a separate guarantee is in existence locally. Pursuant to the following contract terms and conditions, Sigma Dental shall warrant that, at the time of the original purchase, the equipment was free of defects in terms of material, construction and implementation:

1. This limited warranty shall be applicable to the end purchaser of the equipment („Customer“). Any statutory claims of customer or of customer towards seller / dealer shall neither be excluded nor restricted by same.
2. The warranty period shall be twelve (12) months following the purchase of the equipment by the original customer. In case the equipment is resold or otherwise changes owner / user, the warranty shall be valid for the remainder of the 12-month period. These terms and conditions shall not be affected by such change. This limited manufacturer's warranty shall be valid and enforceable only in countries of the European Union.
3. During the warranty period, any defective equipment (at the sole discretion of Sigma Dental) shall be repaired or replaced by Sigma Dental or an authorized customer service. The customer shall receive either the repaired equipment or a fully functional replacement item. Any replaced parts or accessories shall become the property of Sigma Dental.
4. There shall be no extended or renewed warranty period for repaired or replaced equipment.
5. This limited warranty shall not be applicable to any equipment painted or otherwise personalized.
6. Any consumables (such as battery pack) shall have a warranty period of 6 months after first purchase.
7. This limited warranty shall not be applicable to normal symptoms of wear and tear. Furthermore, this limited warranty shall not be applicable if:
 - a) the defect was caused by non-compliance with the operating instructions, improper handling, wetness, moisture, extreme heat or climate conditions or short-term fluctuations of such effects, corrosion, oxidation, unauthorized access or connection attempts, unauthorized opening and repairs, repair attempts with unauthorized spare parts, operator error, improper installation, accidents, natural disasters, spillage of beverages or food, chemical effects or other external effects outside the realm of influence of Sigma Dental (including defects in consumables, such as battery packs, which necessarily have a limited lifespan), unless the defect is directly caused by an error in material, construction, or manufacture;
 - b) the customer does not report the defect to Sigma Dental or the authorized customer service within thirty (30) days after the defect occurred;
 - c) the serial number of the equipment has been removed, scratched out, crossed out or modified or the serial number is illegible;

- d) the defect was caused by the use of the equipment in combination with a part not authorized by Sigma Dental, or by connecting the equipment to such a part, or by using the equipment in some other way in violation of its intended use; (exclusion: USB charging port 5 V)
 - e) the defect was caused by short-circuiting the battery pack, or by breaking the wrapping of the battery pack, or manipulation same, or by using the battery pack for a device other than the intended device.
8. This limited warranty shall represent the only and exclusive legal remedy of the customer against Sigma Dental, and the only and exclusive liability of Sigma Dental toward the customer for defects or functional problems of the equipment. This limited warranty shall replace all other guarantees and warranty statements, whether provided orally or in writing, and whether based on (not necessarily forcefully applicable) statutory regulations, contract provisions or unauthorized acts. Under no circumstances Sigma Dental shall be liable for any ancillary, subsequent or direct damages, costs or expenditures. If the customer is a legal entity, Sigma Dental shall not be liable for any indirect damages, costs or expenditures.
9. Upon assertion of any warranty claim, the original receipt, including the date of sale, must be enclosed.
10. Any changes to this limited warranty shall require the prior written approval by Sigma Dental.

SIGMA DENTAL SYSTEMS-EMASDI GMBH

Heideland 22 | D-24976 Handewitt

Telefon +49 461 - 9 57 88 14

Lichtsysteme@sigmadental.de